


# Microbiological Incubators


**MHP-80** 



Outstanding sample safety for pharmaceutical and clinical laboratories







## Size Info

Width	575mm
Height	760mm
Depth	780mm
Door depth	1086mm
Door width	987mm
Stacking height	80L


Anhui Zhongke Duling Commercial Appliance Co.,Ltd No.1 Hexianghuroad,changfeng,hefei, Anhui, china





### Feature

METHER incubators equip with gradient temperature regulation function and advanced PID algorithm to precisely control the temperature inside the chamber, providing an ideal environment for the cultivation of various microorganisms.






## **Precise Control**



#### Direct Heat & Air-jacketed——

- Rapid temperature recovery in the chamberand reduce ambient temperature interference.
- 4 temperature control zonesare intelligently controlled by microprocessors to ensure temperature uniformity and minimumfluctuation in the chamber .
- Adjustable Temp Range: RT+5~80°C

#### Precise control of temperature Uniformity & Fluctuation



Test points uniformity<±0.3℃</li>

Controlled within ±0.1°C